Fiche Technique

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Références

1. Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
2. Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GIT, 2010.
3. Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

Bibliographie (non exhaustive)

- 1-Mensel B, Wenzel U, Roser M, Lüdemann J, Nauck M. "Considerably Reduced Centrifugation Time Without Increased Hemolysis: Evaluation of the New BD Vacutainer® SST™ II Advance". Clin Chem. 2007; 53(4): 794-5.

2. Smets E and Dijkstra-Lagemaat J. "Influence of Blood Collection in Plastic vs Glass Evacuated Serum-Separator Tubes on Hormone and Tumour Marker Levels". Clin Chem Lab Med, 2004; 42(4): 435-9.
3. Spiritus T, Zaman Z and Desmet W. "Iodinated Contrast Media Interfere with Gel Barrier Formation in Plasma and Serum Separator Tubes". Clin Chem. 2003; 49(7):1187-9.
4. Bakker J, Hackeng C, Church S, Dieijen-Visser M and Beckers O. "An Evaluation of the Integrity of BD Vacutainer® SST™ II and Analyte Stability when Subject to Freezing at -20° C". EUREGIO Congress of Clinical Chemistry and Laboratory Medicine: Aachen, Germany, Oct 2003.
5. Bush V, Janu M, Bathur F, Wells A and Dasgupta A. "Comparison of BD Vacutainer® SST™ Plus Tubes with BD SST™ II Plus Tubes for Common Analytes". Clinica Chimica Acta. April 2001; 306(1-2): 139-43.
6. Bush V, Blennerhasset J, Well A and Dasgupta A. "Stability of Therapeutic Drugs in Serum Collected in Vacutainer Serum Separator Tubes Containing a New Gel (SST II)". Ther Drug Monit. June 2001; 23(3).
7. BD White Paper VS7278-OUS: "A Comparison of BD Vacutainer® SST™ II Advance Tubes with BD Vacutainer® Serum Glass Tubes for Six Infectious Disease Markers". 2006.
8. BD White Paper VS7351: "A Comparison of Adjusted BD Vacutainer® SST™ II Advance Tubes with BD Vacutainer® Serum Glass Tubes for Cortisol, Total T3, Total T4 and TSH on the DPC Immulite® 1000 Analyzer". 2005.
9. BD White Paper VS7050: "Therapeutic Drug Compatibility in BD Vacutainer® SST™ II Plus Tubes". 2004.
10. BD White Paper VS7249: "A Comparative Evaluation of BD Vacutainer® SST™ II Advance Tubes with BD Vacutainer® SST™ Glass Tubes for Select Cardiac Markers". 2004.
11. BD White Paper VS7228: "Performance of BD Vacutainer® SST™ II Advance Tubes at Four and Five Minute Centrifugation Times". 2004.
12. BD White Paper VS7051: "Performance of BD Vacutainer® SST™ II Plus Tubes for Special Chemistry Testing". 2004.
13. BD White Paper VS5780: "Comparison of BD Vacutainer™ SST™ Plus Tubes with SST™ II Plus Tubes Common Analytes on the Toshiba/Abbott Aeroset". 2001.
14. BD White Paper VS5824: "Gel Barrier Stability Comparison of BD Vacutainer® SST™ II Plus, SST™ Plus and SST™ Plus Transport Tubes in Post Centrifugation Transport". 2001.
15. Guder WG, Narayanan S, Wisser H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (Fourth Edition). Darmstadt, Germany: Wiley-VCH; 2009.

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RCS Grenoble B 056 501 711 - Siret 056 501 711 00115. © 2021 BD. BD, le logo BD et toutes les autres marques sont la propriété de Becton, Dickinson and Company.
Becton Dickinson France S.A.S. Société par Actions Simplifiée au capital de 63 356 160 euros.