


Intermittent Female Standard Catheter

Product Data Sheet

Product Codes	See table over page
Description	Standard (Female) Catheters for intermittent bladder drainage.
Intended Use	A single use, flexible tube inserted through the urethra into the bladder for the intermittent/transient passive evacuation of urine, which is then immediately removed. Intended to be used in all patient populations who require bladder drainage. The device will be used by clinicians, trained home carers and patients, in hospitals, clinic, homecare, or any areas for personal use.
Duration of Use	Intermittent but less than 60 minute duration.
Product Information	Standard 23cm in length Short length 18cm Soft rounded tip for atraumatic insertion Two side eyes for maximum drainage
Materials	Tubing and Connector - PVC (DEHP Free)
Quality System	ISO 13485
EU Device Classification	Class Is
Box Qty	See table over page
Packaging	Sterile, individual peel pouches
Sterilised by	Ethylene Oxide;
Shelf Life	5 years
Shore Hardness	80A
Storage & Safety	Toxicity: Manufactured from biocompatible materials. Storage: Store in a clean, cool, dry place away from direct sunlight



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Intermittent Female Standard Catheter

Product Data Sheet

Single Wrapped Product Code	Qty Per Box	Single wrapped Boxed 100	Drug Tariff Code	Qty Per Box	CH/Diam
FC-1406	1000	FC-1406/1	-	-	6ch ø 2.00mm
FC-1408	1000	FC-1408/1	-	-	8ch ø 2.66mm
FC-1410	1000	FC-1410/1	FC-1410/FP	50	10ch ø 3.33mm
FC-1412	1000	FC-1412/1	FC-1412/FP	50	12ch ø 4.00mm
FC-1414	1000	FC-1414/1	FC-1414/FP	50	14ch ø 4.66mm
FC-1416	1000	FC-1416/1	-	-	16ch ø 5.33mm
FC-1418	1000	FC1418/1	-	-	-

Product Code	Qty Per Box	Single wrapped (18cm)	Qty Per Box	CH/Diam
-	-	FC-1406/18/1	100	6ch ø 2.00mm
FC-1408/TRE	450	FC-1408/18/1	100	8ch ø 2.66mm
FC-1410/TRE	450	FC-1410/18/1	100	10ch ø 3.33mm
FC-1412/TRE	450	FC-1412/18/1	100	12ch ø 4.00mm
-	-	FC-1414/18/1	100	14ch ø 4.66mm
-	-	FC-1416/18/1	100	16ch ø 5.33mm
-	-	FC-1418/18/1	100	-

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.