Poole Suckers Suction Devices

Pennine Product Data Sheet

Reference(s)	[See Table]	
Description	Flexible or Rigid sump sucker with removable inner suction probe. Specifically for abdominal surgery, reducing trauma to surrounding tissue.	
Intended Use	Surgical Suction	
Duration of Use	Normally intended for continuous use of less than 60 minutes. Within the time period the duration of use will be patient specific.	
Product Information	Available with multiple side eyes to reduce tissue grab and for quick drainage and suctioning.	
	Removable Inner probe	
	Available with angled tip	
Content Description	Poole Suckers	
Main Materials	YS-3006 & YS-3008 = Flexible PVC (DEHP Free)	
	YS-3007 = SAN	
Quality System	ISO 13485	
EU Classification	Class IIa	
Pack Factor	[See Table]	
Packaging	Sterile, double wrapped individual peel pouches	
Sterilised by	Ethylene Oxide	
Shelf Life	5 years	

Product	Туре	Length	CH/Tip	Inside Diam.	Qty. Per
Code				mm	Box
YS-3006	Flexible straight tip	300	20CH	6.66	50
YS-3008	Flexible (angled tip)	300	20CH	6.66	50
YS-3007	Rigid (White)	300	15CH	5.0	50
-					

STERILE

EO



Storage and Safety:

Toxicity	Manufactured from biocompatible materials.	
Storage	Do not store below 3 °C or above 30 °C	

© 2016 Pennine Healthcare or one of its subsidiaries. The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Revision No 002. Effective Date: Aug 2018





City Gate, London Road, Derby, DE24 8WY, United Kingdom T: +44 (0)1332 794880 www.penninehealthcare.co.uk