Suction Connecting Tubing Surgical Suction


Pennine 📉

Product Data Sheet

Reference(s)	[See Overleaf]			
Description	A complete system to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices. Pennine offer a range of easily connected, interchangeable products in a choice of tubing length and diameter			
Intended Use	The channelling of blood, body fluids and tissues for disposal during surgical suction procedures.			
Duration of Use	Normally intended for continuous use of less than 60 minutes. Within the time period the duration of use will be patient specific.			
Product Information	'F' (Female) connector fits suction machines and male tapered connections for instrumentation			
	'M' (Male) connector push fits into the 'F' connector and is normally used to connect to the standard funnel fitted on suction catheters. Can also be used to connect lengths of tubing together for longer lengths.			
	Application of 'M' connector with 'VC' (Vacuum Control) vents converts a standard catheter to vacuum control application.			
Main Materials	Tube - PVC (DEHP Free)			
	Funnel Connector - PVC (DEHP Free)			
Quality System	ISO 13485			
EU Classification	EU Medical Device Directive Classification I			
Pack Factor	[See Overleaf]			
Packaging	Sterile, individual peel pouches			
Sterilised by	Ethylene Oxide			
Shelf Life	5 years			

STERILE

EO



Shore Hardness and Product CH

Dimensions	[See Overleaf]
Shore Hardness	77
Product CH	[See Overleaf]

Storage and Safety:				
Toxicity	Not applicable, no body contact			
Storage	Store in a clean, cool, dry place, away from direct sunlight.			

© 2016 Pennine Healthcare or one of its subsidiaries. The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



CE

City Gate, London Road, Derby, DE24 8WY, United Kingdom T: +44 (0)1332 794880 www.penninehealthcare.co.uk

Suction Connecting Tube Surgical Suction

Pennine 📉 **Product Data Sheet**

Product Code	Туре	Length	СН	Inside Diameter. mm	Qty Per Box
CT-4312	F/F	1.8m	24CH	5mm	25
CT-4362	F/F	2.7m	24CH	5mm	20
CT-4372	F/F	3.5m	24CH	5mm	20
CT-4314	F/F/M	1.8m	24CH	5mm	25
CT-4364	F/F/M	2.7m	24CH	5mm	20
CT-4334	F/F/M	3.0m	24CH	5mm	20
CT-4394/VC	F/F/M/VC	2.7m	24CH	5mm	20
CT-4374/VC	F/F/M/VC	3.5m	24CH	5mm	20
CT-4314/VC	F/F/M/VC	1.8m	24CH	5mm	25
CT-4324/VC	F/F/M/VC	2.0m	24CH	5mm	25
CT-4334/VC	F/F/M/VC	3.0m	24CH	5mm	20
CT-4012	F/F	1.8m	26CH	6mm	25
CT-4022	F/F	2.0m	26CH	6mm	25
CT-4042	F/F	2.7m	26CH	6mm	20
CT-4032	F/F	3.0m	26CH	6mm	20
CT-4072	F/F	3.5m	26CH	6mm	20
CT4092	F/F	4.0m	26CH	6mm	15
CT-4062	F/F	6.0m	26CH	6mm	10
CT-4020	F/open end	2.0m	26CH	6mm	25
CT-4030	F/open end	3.0m	26CH	6mm	20
CT-4014	F/F/M	1.8m	26CH	6mm	25
CT-4024	F/F/M	2.0m	26CH	6mm	25
CT-4044	F/F/M	2.7m	26CH	6mm	20
CT-4034	F/F/M	3.0m	26CH	6mm	20
CT-4074	F/F/M	3.5m	26CH	6mm	20
CT-4064	F/F/M	6.0m	26CH	6mm	10
CT-4014/VC	F/F/M/VC	1.8m	26CH	6mm	25
CT-4024/VC	F/F/M/VC	2.0m	26CH	6mm	25
CT-4034/VC	F/F/M/VC	3.0m	26CH	6mm	20
CT-4622	F/F	2.0m	30CH	8mm	10
CT-4632	F/F	3.0m	30CH	8mm	10
CT-4672	F/F	3.5m	30CH	8mm	10
CT-4642	F/F	4.0m	30CH	8mm	10
CT-4624	F/F/M	2.0m	30CH	8mm	10
CT-4634	F/F/M	3.0m	30CH	8mm	10
CT-4674	F/F/M	3.5m	30CH	8mm	10
CT-4634/VC	F/F/M/VC	3.0m	30CH	8mm	10

STERILE

EO

Female F/F Female/Female Connectors F/F/M Female/Female/Male Connectors F/F/M/VC Female/Female/Male Connectors with Vacuum Control

F





© 2016 Pennine Healthcare or one of its subsidiaries. The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



CE

City Gate, London Road, Derby, DE24 8WY, United Kingdom T: +44 (0)1332 794880 www.penninehealthcare.co.uk