Oxygen Administration Tubing

Product Data Sheet

Product Codes	See table below		
Description	Oxygen Administration Tubing for connecting oxygen source to patient.		
Intended Use	'Oxygen administration tubes are non-sterile, single use devices intended to provide an additional length of tubing between the oxygen source and patient oxygen delivery device, in patients requiring supplementary oxygen from adults to paediatrics in both clinical, ambulatory, and home care environments. The device will be used by clinicians, trained home carers and patients, in hospitals, clinic, homecare, or any areas for personal use.		
Product Information	Available in non-sterile Coils and Cut Lengths		
	Coils - available in 30m and 50m lengths and supplied in convenient dispenser bag. Funnel ends can be easily cut to fit most oxygen therapy devices.		
	Pre-cut lengths - available in 1.8m and 2.7m lengths, supplied in polythene pack.		
Materials	Tubing - PVC (DEHP Free)		
	Latex Free		
Quality System	ISO 13485		
EU Device Classification	Class IIa		
Box Qty	See table below		
Packaging	Packed non-sterile into individual polythene packs		
Shelf Life	5 years		
Storage & Safety	Toxicity: Manufactured from biocompatible materials		
	Storage: Store in a clean, cool, dry place away from direct sunlight		



Oxygen Administration Tubing - Coil



Oxygen Administration Tubing - Cut Length

Pre-Cut Lengths

Product Code	Length	Qty per box
OC-7311	1.8m	100
OC-7341	2.7m	100

Coil in Dispenser Bags

Product Code	Length	Qty per box
OL-8400	30m	15 Coil
OL-8450	50m	12 Coil

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



