

# Intermittent Standard Nelaton Catheter

## Product Data Sheet

<b>Product Codes</b>	See table over page
<b>Description</b>	Standard Catheters (Nelaton) for intermittent bladder drainage.
<b>Intended Use</b>	A single use, flexible tube inserted through the urethra into the bladder for the intermittent/transient passive evacuation of urine, which is then immediately removed. Intended to be used in all patient populations who require bladder drainage. The device will be used by clinicians, trained home carers and patients, in hospitals, clinic, homecare, or any areas for personal use.
<b>Duration of Use</b>	Intermittent but less than 60 minute duration.
<b>Product Information</b>	38cm in length Soft rounded tip for atraumatic insertion Two side eyes for maximum drainage Available single (SW) or Double (DW) wrapped
<b>Materials</b>	Tubing and Connector - PVC (DEHP Free)
<b>Quality System</b>	ISO 13485
<b>EU Device Classification</b>	Class Is
<b>Box Qty</b>	See table over page
<b>Packaging</b>	Sterile, individual peel pouches
<b>Sterilised by</b>	Ethylene Oxide;
<b>Shelf Life</b>	5 years
<b>Shore Hardness</b>	80A
<b>Storage &amp; Safety</b>	Toxicity: Manufactured from biocompatible materials. Storage: Store in a clean, cool, dry place away from direct sunlight



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

# Intermittent Standard Nelaton Catheter

## Product Data Sheet

Single Wrapped Product Code	Qty Per Box	Drug Tariff Code	Qty Per Box	Single Wrapped Boxed 100	CH/Diam
NC-1206/SW	500	NC-1206/FP/25	200	NC-1206/SW/1	6ch ø 2.00mm
NC-1208/SW	500	NC-1208/FP/25	200	NC-1208/SW/1	8ch ø 2.66mm
NC-1210/SW	500	NC-1210/FP/25	200	NC-1210/SW/1	10ch ø 3.33mm
NC-1212/SW	500	NC-1212/FP	200	NC-1212/SW/1	12ch ø 4.00mm
NC-1214/SW	500	NC-1214/FP	200	NC-1214/SW/1	14ch ø 4.66mm
NC-1216/SW	500	NC-1216/FP	200	NC-1216/SW/1	16ch ø 5.33mm
NC-1218/SW	500	-	-	NC-1218/SW/1	18ch ø 6.00mm
NC-1220/SW	400	-	-	NC-1220/SW/1	20ch ø 6.66mm
NC-1222/SW	400	-	-	NC-1222/SW/1	22ch ø 7.33mm
NC-1224/SW	400	-	-	-	24ch ø 8.00mm

Double Wrapped Product Code	Qty Per Box	Double Wrapped Inner bag 25s	Qty Per Box	CH/Diam
NC-1210/DW	400	NC-1210/DW/TRE	400	10ch ø 3.33mm
NC-1212/DW	400	NC-1212/DW/TRE	400	12ch ø 4.00mm
NC-1214/DW	400	NC-1214/DW/TRE	400	14ch ø 4.66mm
NC-1216/DW	400	-	-	16ch ø 5.33mm

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.