

Reference(s)	[See Overleaf]
Description	Rigid, transparent suction devices for use with suction connecting tubing. For general purpose surgical suctioning. Also available as a set, pre-attached to 2 or 3 metre tubing.
Intended Use	Surgical Suction
Duration of Use	Normally intended for continuous use of less than 60 minutes. Within the time period the duration of use will be patient specific.
Product Information	Available with or without vacuum control (VC) vent Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip Midi Yankauer available with option of flat tip or Atrauma Tip to reduce tissue grab.
Content Description	Link Yankauer
Main Materials	PVC (DEHP Free)
Quality System	ISO 13485
EU Classification	Class IIa
Pack Factor	[See Overleaf]
Packaging	Sterile, double wrapped individual peel pouches
Sterilised by	Ethylene Oxide
Shelf Life	5 years



Shore Hardness and Product CH

Dimensions [See Overleaf]

Product CH [See Overleaf]

Storage and Safety:

Toxicity Manufactured from biocompatible materials.



Storage Do not store below 3 °C or above 30 °C



© 2016 Pennine Healthcare or one of its subsidiaries.

The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Product Code	Type	Use	Length	CH	Qty. Per Box
LY-3400	Short Mini	Paediatric/dental (with 2 side eyes)	190mm	15ch	100
LY-3401	Short Mini VC		190mm	15ch	100
LY-4400	Short Mini Supertip		195mm	15ch	100
LY-4401	Short Mini Supertip VC		195mm	15ch	100
LY-4420	Short Mini Super Tip		195mm	15ch	100
LY-3500	Mini Fine Tip	Fine suctioning (e.g. ENT)	270mm	12ch	50
LY-3501	Mini Fine Tip VC		270mm	12ch	50
LY-4500	Mini Fine Supertip		275mm	12ch	50
LY-4501	Mini Fine Supertip VC		275mm	12ch	50
LY-3600	Midi 30° angle	For general and pharyngeal suction	280mm	18ch	50
LY-3601	Midi VC		280mm	18ch	50
LY-3600/PR	Midi 10° angle		270mm	18ch	50
LY-4600	Midi Supertip		280mm	18ch	50
LY-4601	Midi Supertip VC		280mm	18ch	50
LY-3700	Maxi	For large particles and significant volumes of fluid	275mm	22ch	50
LY-3701	Maxi VC		275mm	22ch	50
LY-4700	Maxi Supertip		275mm	22ch	50
LY-4701	Maxi Supertip VC		275mm	22ch	50
LY-3800	Maxi Plus	For orthopaedic use	280mm	28ch	50
LY-3801	Maxi Plus VC		280mm	28ch	50
LY-3602	Atrauma tip, short soft	For general suction with minimal tissue grab	280mm	18ch	50
LY-3603	As above with VC		280mm	18ch	50
LY-3606	Atrauma tip, short rigid	For general suction with minimal tissue grab	280mm	18ch	50
LY-3607	As above with VC		280mm	18ch	50



© 2016 Pennine Healthcare or one of its subsidiaries.

The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.