# Flexi-Rib Tubing Suction Connecting Tubing



### **Product Data Sheet**

Reference(s)	[See Overleaf]				
Description	Flexi-Rib tubing to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices.				
Intended Use	The channelling of blood, body fluids and tissues for disposal during surgical suction procedures.				
Duration of Use	Normally intended for continuous use of less than 60 minutes. Within the time period the duration of use will be patient specific.				
	Ribbed connecting tubing for enhanced drapability				
	Non-slip outer rib reduces friction				
Product Information	Resistant to coiling, kinking				
	Maximum pressure 500mm HG Mercury				
	Colour coded funnel for diameter identification				
	Multi-fit to all suction devices and sources				
	Economical box size				
Main Materials	Tube - PVC (DEHP Free)				
- Walli Materials	Funnel/Connector - PVC (DEHP Free)				
Quality System	ISO 13485				
EU Classification	EU Medical Device Directive Classification I				
Pack Factor	50 per box				
Packaging	Sterile, double wrapped individual peel pouches				
Sterilised by	Ethylene Oxide				
Shelf Life	5 years				



#### Shore Hardness and Product CH

Dimensions	[See Overleaf]		
Shore Hardness	71A		
Product CH	[See Overleaf]		

### Storage and Safety:

Toxicity	Manufactured from biocompatible materials
Storage	Store in a clean, cool, dry place, away from direct sunlight.









© 2016 Pennine Healthcare or one of its subsidiaries.

The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



# Flexi-Rib Tubing Suction Connection Tubing



## **Product Data Sheet**


Product Code	Туре	Length	CH	Inside Diameter. mm	Qty Per Box
CT-4922	F/F	2.0m	26CH	6	60
CT-4932	F/F	3.0m	26CH	6	50
CT-4972	F/F	3.5m	26CH	6	50
CT-4924	F/F/M	2.0m	26CH	6	60
CT-4934	F/F/M	3.0m	26CH	6	50
CT-4722	F/F	2.0m	28CH	7	60
CT-4732	F.F	3.0m	28Ch	7	50
CT-4772	F/F	3.5m	28CH	7	40
CT-4724	F/F/M	2.0m	28CH	7	60
CT-4734	F/F/M	2.0m	28CH	7	50
CT-4774	F/F/M	3.5m	28CH	7	50
CT-4714/VC	F/F/M/VC	1.8m	28CH	7	25
CT-4741/VC	F/F/M/VC	2.5m	28CH	7	20











The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

