5-in-1 Connector Suction Connecting Tubing



Product Data Sheet

Reference(s)	[See Table]	
Description	A 5-in-1 connector to connect two tubes together.	
Intended Use	To connect surgical tubing of different diameter	
Duration of Use	Normally intended for continuous use of less than 60 minutes. Within the time period the duration of use will be patient specific.	
Product Information	Can be used to extend lengthof tubing (i.e to connect two tubes together of either the same or different diameters)	
	Designed for use with all types of tubing	
	Available either bulk, non sterile or packed sterile	
Content Description	5-in-1 Connector	
Main Materials	PVC (DEHP Free)	
Manufacturer	Pennine Healthcare ISO 13485	
Regulation	EU Medical Device Directive Classification STERILE	
Pack Factor	[See Table]	
Packaging	Sterile, individual peel pouches	
Sterilised by	Ethylene Oxide	
Shelf Life	5 years	

		Qty. Per Box
DP-501	Sterile	100
DP-501/Bulk	Bulk, non Sterile	500



Toxicity	Manufactured from biocompatible materials
Storage	Do not store below 3 °C or above 30 °C









© 2016 Pennine Healthcare or one of its subsidiaries.

The information given in this sheet is not intended to be exhaustive. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

