


Product Specification



Article	OR Gown, <i>Mediwell light</i> <i>Standard Performance</i>
Specification No.	46 44 33C
Article No.	46 44 33 <i>sterile goods</i> 46 44 33/9 o.E. , <i>non-sterile goods</i> 46 44 33/9 m.E. , <i>non-sterile goods</i>
1. Description	OR Gown made of Mediwell light , <i>non reinforced</i> . The OR Gown and two hand towels are supplied in an overwrap in <i>sterilisation packaging</i> .
2. Application	For procedures without a particularly high fluid incidence.
3. Production	The above mentioned product is manufactured in hygienic conditions.
4. Material composition	<i>all of the materials used are <u>latex-free!</u></i>
4.1 <u>Mediwell light</u>	
• Nonwoven	polypropylene nonwoven, turquoise
4.2 Cuff	synthetic fibre
4.3 <i>Included loose</i>	
• Hand towel 2x	# 62 40 40
5. Product properties	fluid and alcohol repellent, high air permeability, highly comfortable to wear, low-lint.
6. Technical data	
6.1 <u>Dimensions</u>	
• Size	large
• Length	130 cm ± 5 cm
6.1.1 <i>Cuff</i>	
• Length	8 cm ± 2 cm
6.1.2 <i>Finale folded dimensions</i>	17 x 20 cm ± 3 cm
6.1.3 <i>Included loose / Hand towel 2x</i>	
• Length	39.5 cm ± 2 cm
• Width	30 cm ± 1 cm

Quality Assurance:	
Date: 20.09.2022	Signature: 



Article OR Gown, *Mediwell light*
Standard Performance

Specification No. 46 44 33C

6.2 Raw material

- *Nonwoven*
- Basis weight 35 g/m² ± 10 %

6.3 Performance requirements

- The material used to make this product exceeds the minimum requirements, in accordance with the pertinent parts of EN 13795.
- Individual values can be found in the corresponding material specification.

The raw materials used in manufacture are physiologically harmless.


Product Specification



Article OR Gown, *Mediwell light*
Standard Performance

Specification No. 46 44 33C

7. Schematic diagram




all dimensions in cm



Article OR Gown, *Mediwell light*
Standard Performance

Specification No. 46 44 33C

7.1 Folding




all dimensions in cm



Article OR Gown, *Mediwell light*,
Standard Performance

Specification No. 46 44 33C

7.2 Overwrap



all dimensions in cm



Article	OR Gown, <i>Mediwell light Standard Performance</i>
Specification No.	46 44 33C
8. Packaging / sterile goods	
8.1 <i>Primary packaging</i>	
• Paper	60 g/m ² ± 10 %
• Film	100 µm ± 10 %, transparent
8.2 <i>Film bag</i>	
• Size	990 x 800 mm ± 10 mm
• Film	PE film, transparent
• Contents	32 pcs.
8.3 <i>Transport carton</i>	
• Size	600 x 400 x 260 mm ± 10 mm
• Contents	32 pcs.
8.4 Sterilisation	EtO sterilisation (<i>in accordance with ISO 11135</i>)
8.5 <i>Usability and storage</i>	5 years, no direct exposure to sunshine, store at a dry place
9. Packaging / non-sterile goods	
9.1 <i>Film bag 2x</i>	
• Size	1.200 x 800 mm ± 10 mm
• Film	PE film 35 µm ± 10 %, transparent
9.2 <i>Transport carton</i>	
• Size	780 x 385 x 295 mm ± 10 mm
• Contents / # 464433-9o.E.	60 pcs./ <i>without overwrap</i>
• Contents / # 464433-9m.E.	50 pcs./ <i>with overwrap</i>
9.3 Sterilisation	suitable for sterilisation by EtO gas.
9.4 <i>Usability and storage</i>	6 years, no direct exposure to sunshine store at a dry place
10. Disposal	The raw material can be disposed of in an environmentally friendly manner. No toxic gases are released during incineration. The ground water is not contaminated by toxic substances when dumped.

This product specification is the property of SENGEWALD Klinikprodukte GmbH and is to be treated with the strictest confidence. Communication to third parties and duplication, in whole or in part, are not permitted.