



Medi-Vac Guardian Fluid Management Suction Canisters

Technical Sheet

Reference (s)	SKU	Long description	Volume	UOM
	65651-212	Medi-Vac rigid suction canister liner with mechanical filter	1200 mL	100 ea/cs
	65651-220	Medi-Vac rigid suction canister liner with mechanical filter	2000 mL	40 ea/cs
	65651-230	Medi-Vac rigid suction canister liner with mechanical filter	3000 mL	40 ea/cs
	65651-412	Medi-Vac rigid suction canister liner with PPV filter	1200 mL	40 ea/cs
	65651-420	Medi-Vac rigid suction canister liner with PPV filter	2000 mL	40 ea/cs
	65651-430	Medi-Vac rigid suction canister liner with PPV filter	3000 mL	40 ea/cs

Description	Suction canisters are used to receive the body fluid removed from the patient. The collected liquids are not intended to be reinfused into the body nor reprocessed for eventual reinfusion, but simply discarded.
Main Materials	Lid: Polypropylene Canister: Polystyrene
Manufacturing	Address: Cardinal Health 1430 Waukegan Rd Waukegan, IL 60085 Phone: +1.847.473.1500
Regulation	CE marked. Single use only. Consult instructions for use. 2007/47/EC, EN ISO 14971:2012, ISO 15223-1:2012, EN 1041:2008, 93/42/EEC, 80/181/EEC, ASTM D4169-14, ISO 5359:2014, ISO 13485:2012, ISO 11607-1:2006, ISO 10993-1:2009, ISO 30993, EN 10079-3:2014
Sterilization type	Canisters and accessories are offered as non-sterile devices only

cardinalhealth.com