


Product category: Sterile Syringes with NaCl 0,9%

Description: Syringe NaCl 0,9% Luer Lock 10ml Sterile/Non Sterile

Technical drawing:



Material: Syringe made of polypropylene
Syringe cap made polypropylene
Sodium chloride 0,9%

REF	Designation	Dimensions (L x W x H)
14360SNS	Syringe NaCl 0,9% Luer Lock 3 ml Sterile/Non Sterile	240 x 50 x 2mm
14361SNS	Syringe NaCl 0,9% Luer Lock 5 ml (3 ml filling) *S & NS Out	210 x 50 x 2 mm
14362SNS	Syringe NaCl 0,9% Luer Lock 5 ml Sterile/Non Sterile	210 x 50 x 2 mm
14363SNS	Syringe NaCl 0,9% Luer Lock 10 ml Sterile/Non Sterile	210 x 50 x 2 mm
14364SNS	Syringe NaCl 0,9% Luer Lock 20 ml Sterile/Non Sterile	270 x 65 x 25 mm
14365SNS	Syringe luer lock 10 ml NaCl 0,9% (3ml fill) sterile&NS	240 x 50 x 2 mm
14366SNS	Syringe luer lock 10 ml NaCl 0,9% (5ml fill) sterile& NS	240 x 50 x 2 mm

Sterilization: Sterile fluid path by radiation e-beam

Brand: Steriflush

Product overview: Single pack

REF	Description	Quantity in pcs	
		Inner box	Transport carton
14360SNS	Syringe NaCl 0,9% Luer Lock 3 ml Sterile/Non Sterile	150	600
14361SNS	Syringe NaCl 0,9% Luer Lock 5 ml (3 ml filling) *S & NS Out	150	600
14362SNS	Syringe NaCl 0,9% Luer Lock 5 ml Sterile/Non Sterile	150	600
14363SNS	Syringe NaCl 0,9% Luer Lock 10 ml Sterile/Non Sterile	125	500
14364SNS	Syringe NaCl 0,9% Luer Lock 20 ml Sterile/Non Sterile	60	240
14365SNS	Syringe NaCl 0,9% Luer Lock 10 ml (3 ml filling) *S & NS Out	125	500
14366SNS	Syringe NaCl 0,9% Luer Lock 10 ml (5 ml filling) *S & NS Out	125	500



Sterisets Medical Products
Zona Industrial 1 - Lote 11 a 14
4560-164 Guilhufe, Penafiel
Portugal www.sterisets.com

Conforms to ISO 13485, Directive 93/42/EEC
and Regulation (EU) 2017/745

Effective: 18-03-2021

Issued: DQ

Approved: DQ



REF 14360SNS
& variants

MDDTF_L-4_11


ET 1308/9

IMP 145/4

Cap drawing:





剖面 A-A



Label drawing:

6.38mm
3.07mm
Pre-printed transparent label



Sterisets Medical Products
Zona Industrial 1 - Lote 11 a 14
4560-164 Guilhufe, Penafiel
Portugal www.sterisets.com

Conforms to ISO 13485, Directive 93/42/EEC
and Regulation (EU) 2017/745

Effective: 18-03-2021

Issued: DQ




Approved: DQ

REF 14360SNS
& variants

MDDTF_L-4_11

ET 1308/9

IMP 145/4



Pre-printed transparent label

